

The refined choice



The Fine Range - accuracy you demand and quality as standard

Complimenting the traditional Swann-Morton blade shapes and developed potentially for more complex and challenging surgical procedures which may be encountered in such disciplines as reconstruction, ophthalmics, cardiology and even hair restoration, the surgeons may sometimes prefer our Fine Range.

The Fine Range has small, slim, precise blades in stainless steel for delicate operations.

Six special handles, including the SF13 and SF23 with independently screwed top sections to facilitate

blade-changing, take all blades, in stainless steel, sterile, individually sealed in metal foil packets and boxed in units of 25.

Handles are supplied individually in plastic wallets.

Sterile blades should be selected for all invasive surgical procedures and are for single patient use thus minimising the risk of cross infection.

Swann-Morton®
SHEFFIELD ENGLAND



**UK
CA**
All medical devices comply with the Medical Device Directive, 93/42/EEC as referred by Part II of The Medical Devices Regulations 2002, Annex 2, excluding section 4 (as modified by Part 2 of Schedule 2A to The Medical Devices Regulations 2002)



Scan this code for more product information




ONE OF THE WORLD'S LEADING SURGICAL BLADE MANUFACTURERS

Product	Product REF	Unit	Product	Product REF	Unit
Fine Range Blades - Sterile					
Blades in Stainless Steel - Sterile, individually V.P.I. protected in Metal Foil Packets					
Fits handles SF1, SF2, SF3, SF4, SF13 & SF23			Fits handles SF1, SF2, SF3, SF4, SF13 & SF23		
SM61	5901		SM61(SB)	5911	
SM62	5902		SM62(SB)	5912	
SM64	5904		SM63	5903	
SM65	5905	25 Blades individually packed	SM69	5909	25 Blades individually packed
SM65A	5906		SP90	5921	
SM67	5907		SP91	5922	
SM68	5908				
			Myringotomy	2102	10 Blades individually packed

Specialist Applications

Cardio-Vascular Surgery: SM65, SM65a, SM67. Plastic Surgery: SM63, SM64, SM67 & SM69. Ophthalmology: SM63, SM64, SM67 & SM69. Orthopaedic Surgery: SM64, SM65 & SM65a. E.N.T. Surgery: SM68 & Myringotomy. Podiatry: SM61, SM61(SB), SM62, SM62(SB), PD81 & PD82. Dental Periodontics & Implantology: SM65, SM67 & SM69. Hair Restoration & Transplantation: SP90 & SP91

Product	Product REF	Unit	Product	Product REF	Unit
Fine Range Handles					
1 handle individually packed					
1 handle individually packed (2 piece)			1 handle individually packed (3 piece)		
SF1	6051	1	SF13	6061	1
SF2	6052	1	SF23	6062	1
SF3	6053	1			
SF4	6054	1			



Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.