




Safe



their feet in
your hands



Swann-Morton the preferred choice of Podiatrists and Chiropodists.

Swann-Morton has become a world leader in the production of quality surgical blades, scalpels and handles and the preferred choice of Podiatrists and Chiropodists through their precision, consistency, and reliability.

A wide and varied range of products provide the practitioner with a choice of tools for each treatment stream enabling them to optimise their patient outcomes and satisfaction.

Swann-Morton®
SHEFFIELD ENGLAND



All medical devices comply with the Medical Device Directive, 93/42/EEC as referred by Part 2 of The Medical Devices Regulations 2002. Annex A (excluding section 4) as modified by Part 2 of Schedule 2A to the Medical Devices Regulations 2002.



Scan this code for more product information



ONE OF THE WORLD'S LEADING SURGICAL BLADE MANUFACTURERS

Overtime we can develop corns and calluses in high pressure areas of the feet such as the metatarsal heads, the sides and tops of toes and the heel which require regular pairing down by a trained Podiatrist.

Podiatrists along with their colleagues in the NHS are currently faced with meeting the challenges presented by recent legislation on the prevention of sharps injuries amongst Healthcare workers.


Training in the best practice for the safe handling of sharps has become essential and guidelines should be followed whether in a private practice, GP practice or a hospital setting.

Whether in-practice or mobile in the community the new versatile "BladeFlask EVO" embraces safety protocols and is compliant with UK Health and Safety (Sharp Instruments in Healthcare) Regulations 2013. A sealed puncture proof container that can quickly, efficiently, and safely remove a contaminated blade from an extensive range of handles including the hexagonal Barron (B3) to a Bulbous 5 or 6 in a single-handed movement that can be conveniently located within each practice to safely remove and store up to 100.

In addition, the Swann-Morton Blade Remover allows for safe and easy removal of all sizes of scalpel blade from both No.3 and No.4 handles.

Supplied sterile the transparent nature of the material allows for 100% accountability of the blades used during a procedure, each one safely encapsulated within an individual remover which may then be disposed of in an appropriate sharps bin or container.

A Step forward in foot care



The Swann-Morton Sabre E/11 blade remains popular amongst Podiatrists for removing those areas of hard skin thus minimising discomfort. The special "closed" razor like cutting edge is ideal for this procedure. Nail care may be needed for thick painful and difficult to manage toenails, ingrowing nails, and fungal nails whilst Nail Surgery may be the only

permanent solution for painful, ingrowing, involuted or thickened deformed toenails. The Fine SM61 and SM62 "single bevel blades" act a like a small chisel which can be inserted in and around the nail whilst the PD81 blade used in conjunction with the slim barrelled PD handle offer an extra level of "finish" to the nail.

More information can be found on the Swann-Morton website at www.swann-morton.com/podiatry.php

Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.